



4NRJ
WWW.4NRJ.COM

CONNECTED CLINALERT MINI

NEW VERSION



FEATURES

Ref : CLINALERTMINI



- ✓ Network : LoRaWAN®
- ✓ Utilisation : all metal supports
- ✓ Fixation device : magnet (4NRJ patent)
- ✓ Temperatures reading : +/- 0.5 C°
- ✓ Autonomy up to 15 years
- ✓ Accuracy : +/- 0,001°
- ✓ Alimentation : Lithium battery : 3.6V
- ✓ Protection rating : IP67
- ✓ Dimensions : 105 x 79 x 49 mm
- ✓ Mass weight : 0.57 kg

Setting tool, Ref : CLINALERTMINIOP

- ✓ Fixation on catenary pole by P-FIX system (4NRJ patent)

Application 4NRJ

- ✓ CLINALERT Download free



01

CLINALERT MINI solution allows a continuous information of the inclination of the support to watch the verticality overtime. Inclination thresholds can be customized and alerts sent to SMS and/or E-mail in case of exceeding thresholds.




02

CLINALERT MINI solution is set thanks to a dedicated insulating pole and related setting tool. The idea is to position the CLINALERT device at an appropriate height to avoid thefts or damage.

03

Download 4NRJ free application CLINALERT MINI available in Play Store®. A simple contact between your Smartphone and your CLINALERT MINI device is enough to activate the application.




CLINALERT MONITORING AND MEASUREMENTS READING ARE AVAILABLE FROM YOUR PC

- Data recording
- Follow up all support inclinations at any time night and day.
- Periodic monitoring since the installation of the CLINALERT solution.
- Configurable display to analyse the desired timeframe.
- Remotely reset to zero if necessary.
- Remotely change the alerts thresholds.
- The initial position of the post is the reference value.



REMOTELY WATCH
THE VERTICALITY
OF YOUR CATENARY
POSTS



4NRJ
WWW.4NRJ.COM

2 Rue Albert Calmette
ZA Les Gailletrous
41260 La Chaussée-St-Victor
France

Tél : +33 (0)2 54 42 05 12
E-mail : contact@4nrj.com